

Osteotec

Silicone Finger Implant

The Osteotec Silicone Finger Implant is a one piece, flexible silicone elastomer implant placed during arthroplasty into the proximal and distal intramedullary canals of either the MCP or PIP joint.

Over **20 years of clinical use**

Available sterile in eleven sizes (00–9)


Studies demonstrate **strong clinical survivorship at 10 years**, with associated improvements in key patient reported outcomes

Recorded **over 17,000 sales** of Osteotec Silicone Finger Implants since its introduction to the market in 2004

Instrumentation design

A single-tray instrument set enables efficient preparation and sizing:

- Comprised of the following, all supplied non sterile:
 - Starter awl
 - Proximal and distal rasp (sizes 00–9)
 - Osteotec Silicone Finger Implant Sizers (sizes 00–9)



Scan to view our animation, product page and website:

